

Polarization controllers: key differentiators to optimize your testing and manufacturing applications

Sponsored by: Luna Innovations

After optical power, polarization is perhaps the next most important test and measurement optical property of fiber optic components and systems. Specifically, with the introduction of coherent dual-polarization fiber optic communications, fiber optic sensing, integrated optics and quantum communications and computing, polarization control has become an integral part of optical R&D labs and manufacturing lines. In contrast to manual polarization controllers, multifunction polarization controllers (MPCs) employ electronic control of optical polarization actuators to convert any given state of polarization (SOP) to any other desired SOP. The waveplates within the MPC components can be made from various materials and technologies including magneto-optics birefringent crystals, liquid crystals, piezo-electric fiber squeezers or LiNbO₃. Two complimentary MPC technologies covering a wide range of applications are detailed below.

After optical power, polarization is perhaps the next most important test and measurement optical property of fiber optic components and systems.

Fiber squeezers advantages and considerations

There are several technologies available to realize electronic polarization control. Among them, fiber squeezers mechanically compress the fiber creating birefringence. This “squeezing” of the fiber causes one polarization of light to be retarded with respect to its orthogonal polarization. Often, three or four voltage-controlled squeezers are sequentially placed at different orientations along the fiber to enable total control of the SOP, based upon the voltages applied (Figure 1).

Luna’s MPC 201, MPC 202 and MPC 203 combine fiber squeezer polarization controllers with proprietary polarization control algorithms to achieve a wide range of polarization control functions specifically scrambling the SOP and setting the SOP functions.




Figure 1. Fiber squeezing technology physically compresses the fiber with a series of electrically driven piezo actuators placed at alternating orientations around the fiber for complete control.
Source: Luna Innovations.

In addition to the squeezers’ seemingly simple method, this technology has other considerable advantages. Fiber squeezers are all-fiber devices and are as broadband as the optical fiber, making squeezers available for any wavelength window. For instance, the most common squeezers, employing standard SMF fiber, operate across the standard datacom and telecom windows, 1260 nm to 1650 nm. Squeezers also have extremely low insertion loss, low back reflection and low polarization dependent loss (PDL). These characteristics make squeezer MPCs useful for high-precision PDL instrumentation.

Sponsored by:



Produced by:



Depending on the application, squeezers may be not always be the best choice. The mechanical compression of the fiber causes the squeezer technology to have a response time in the μsec to msec range. Furthermore, the birefringence change is limited by the maximum pressure that can be applied to the fiber. As such, the SOP is not endlessly controllable. For applications where speed and endless polarization control are desired, Luna's LiNbO_3 polarization controllers are a better choice and compliment the squeezer product line.

LiNbO₃ controller advantages and considerations

In Luna's LiNbO_3 polarization control products, the optical LiNbO_3 chip device has a defined waveguide to guide the light from the input fiber to the output fiber. A set of three electrodes straddling the waveguide form an electro-optic waveplate. Properly applying voltages to the electrodes will control a variable birefringence endlessly rotatable waveplate. LiNbO_3 controllers can have up to eight waveplates, as shown in Figure 2, each with both birefringence and angular waveplate control, totaling 16 degrees of freedom. Not only does the LiNbO_3 have a fast electro-optic response where the polarization can be modified in less than a microsecond, but the polarization control is endless. These underlying features allow for the implementation of flexible, ultra-fast, compound and complicated polarization control functions. The newly released MPC 2500 provides seven different polarization functions.




Figure 2. Luna's ultrafast polarization control device driven by high slew-rate voltage drive amplifiers and drive circuitry to run tracking algorithms and polarization control processes. Source: GlobalSpec/Luna Innovations

Both of Luna's polarization control technologies have functions to either control the SOP or mix up the SOP. MPC functions for stress and verification testing, include Luna's "Scrambler," "Spinner," "Depolarizer" and "Randomizer," each mixes up the SOP differently for specific testing purposes. These SOP functions serve both general lab uses and telecom system testing. Functions that set the SOP include "Tracker" or "Stabilizer," "Acquirer" and "Paddle" mode drive the SOP to a desired state. These modes serve more specific applications like silicon photonics testing and coherent beam combination. Table 1 includes all functional modes of operation of the MPC 2500 LiNbO_3 .

The MPC 2500 is an updated and rebranded version of the NRT-2500, and joins the MPC 201, MPC 202 and MPC 203, squeezer-based test and measurement products, in a single multifunction polarization controller portfolio.

The LiNbO_3 offers more functions than its fiber squeezer counterparts; however, a deeper consideration of the end-user's application requirements is crucial to properly select the optimal MPC product.

Function		Description	Use Case
Mix up SOP	Scrambler	Stochastic Rayleigh scrambling of SOP uniformly covering the Poincaré sphere with a quantifiable statistical speed distribution.	Stress testing polarization insensitivity or tracking capabilities in lab and manufacturing application
	Depolarizer	An extreme scrambler to move SOP faster than a detector response time, removing the effects of polarization.	Remove PDL from a system
	Spinner	Generates a single speed of SOP change.	Determines polarization demultiplexing capabilities of dual polarization coherent DSP receivers
	Randomizer	A sudden/random jump of SOP	Simulate the SOP temporal response to lightning strikes on coherent DSP receivers
Control the SOP	Paddles	General use manual adjustment to set the SOP at a desired state.	Generic
	Acquirer	A "smart" paddle mode uses feedback to automatically drive SOP to desired state.	Used in Si photonics testing to test polarization dependent loss or polarization dependent coupling
	Tracker or Stabilizer	Continuous manipulation of moving polarization inputs or outputs to stabilize SOP at a single state.	Used for coherent signal and LO alignment and coherent beam combination

Table 1. all functional modes of operation of the MPC 2500 LiNbO_3 , Source: Luna Innovations.

Fiber squeezers and LiNbO₃ controller applications – Two complementary technologies

Table 2 shows a more direct comparison between these technologies; one technology is strong when the other is less capable. Ultimately these tradeoffs will impact which one of these devices is best suited for the end-user’s application.

	MPC 203 (Squeezer)	MPC 2500 (LiNbO ₃)
Insertion loss	< 0.5 dB*	< 3 dB
PDL	< 0.05 dB	< 0.3 dB
Optical power handling	30 dBm*	20 dBm
Bandwidth	1260 nm to 1650 nm	requires more calibration
Device type	all fiber	fiber-pigtailed waveguide
Birefringent polarization control	yes	yes
Angular polarization control	no	endless rotation
Phase control	no	yes
Tracking response time	10-3 seconds	10-6 seconds

Table 2. A direct comparison between MPC 203 and MPC 2500. *Dominated by connectors. Source: Luna Innovations

General lab use versus dynamic environments

The MPC response time must also be considered when deciding which polarization control technology to pick. SOP fluctuations can occur in a lab, or naturally occurring from wind and weather (e.g. aerial fiber networks), trains, cars, trucks and heavy equipment passing fibers in close proximity to fiber installation (e.g. buried fiber networks). While a single fiber squeezer can change the SOP in about 30 μs, tracking and stabilization require multiple SOP steps to find and lock onto the desired SOP. Mechanical squeezers, such as the MPC 201, MPC 202 and MPC 203, are generally fast enough to track or mitigate most typical mechanically induced SOP fluctuations in the lab up to a few radians per second. LiNbO₃ controllers are usually preferred when dynamic SOP fluctuations are faster than a few milliseconds, due to their significantly faster response time.

Mechanical squeezers, such as the MPC 201, MPC 202 and MPC 203, are generally fast enough to track or mitigate most typical mechanically induced SOP fluctuations in the lab up to a few radians per second. LiNbO₃ controllers are usually preferred when dynamic SOP fluctuations are faster than a few milliseconds, due to their significantly faster response time.

Multi-band operation

Polarization scrambling and polarization setting functions are often required for fiber optic component testing across the optical windows of single mode fiber, especially in the O-band (1260-1360 nm) and L-band (1530 – 1565 nm). Both the LiNbO₃ (MPC 2500) and the squeezer (MPC 20x) MPCs work well, and can be used to scan through all possible polarization for adequate testing and characterization. The exception is the Spinner function of the MPC 2500, which needs more precise wavelength calibration of the waveplates. The intrinsic wide wavelength coverage of the all-fiber squeezer MPCs can simplify test station and reduce cost.

Quantum

Quantum communications applications, such as quantum key distribution, take advantage of the laws of quantum physics to protect data. In some implementations, the quantum channel uses a single photon in a superposition polarization state (qubit) to identify a hacker’s presence on the optical fiber or free space communication channel. Quantum computing employs quantum entangled qubits to enable calculations intractable by normal (i.e. classical) computers. MPCs are often used here to prepare the entangled qubits. In both quantum communications and computing, the optical quantum channel has only a few photons to carry the information. Hence loss in the fiber channel is intolerable. For these applications the low loss squeezer is the preferred MPC technology.

Coherent beam combination

Coherent beam combination is another relatively nascent field where amplified clones of a single laser are re-combined to create high-power directed energy laser sources for military applications, welding or cutting of metal parts. A number of techniques can be employed to scale up the power to maximal level. These techniques generally favor the ultra-fast response of LiNbO₃ controllers to combine the substantial number of beam elements. The Tracker function on the MPC 2500 is optimized for these applications where the active tracking algorithm allows for a nearly perfect SOP and phase tracking, enabling the system to coherently align all beams to the master beam.

Silicon photonics

As shared in Table 1, silicon photonics testing and telecom applications will require tests where polarization dependent loss or polarization dependent coupling needs to be evaluated. In these cases users may desire to rapidly drive the SOP to maximize or minimize a feedback signal. The ultra-fast LiNbO₃ polarization control – driven by high slew rate voltage drive amplifiers – quickly responds to optical feedback and directs the change in SOP in milliseconds with the Acquirer function. The high precision and fast response of the MPC 2500 is preferable to increase the throughput of optical wafer and die tests in a Si photonic chip fab.

Telecom

State-of-the-art, high bit rate communications systems, with coherent polarization-multiplexed modulations, must undergo receiver verification testing; this ensures the polarization demultiplexing capabilities are sufficient to track changes of the network SOP. The Spinner function of the MPC 2500 generates a well-defined speed of change for the SOP, rotating an electrooptic ½-waveplate up to 1 Mrad/sec (Figure 4). A second application function that exploits the high speed of the MPC 2500 is the Randomizer. It randomly jumps the SOP in 10 – 20 μs to emulate the ultrafast SOP changes caused by lightning along aerial fiber optic links.




Figure 3. Sample histogram of spinner function. Source: News Ridge Tech

Conclusion

MPCs can transform any arbitrary input SOP to any output SOP and create dynamic SOP functions to test, stress and evaluate optical devices. The versatility of these devices cannot be understated. However, the different technologies used in MPCs can potentially lead to some confusion when it comes to deciding what type of MPC to employ. It may seem as though fiber squeezer MPCs are in direct competition with LiNbO₃ technologies; this is not necessarily the case, as the advantages and considerations for each complement each other more than they conflict with each other.

	MPC 203 (Squeezer)	MPC 2500 (LiNbO ₃)
General lab use	recommended	
Broadband scanning	recommended	
Quantum applications	recommended	
Dynamic environments		recommended
Coherent beam combination		recommended
Silicon photonics		recommended
Telecom		recommended

Table 2. A direct comparison between MPC 203 and MPC 2500.
*Dominated by connectors. Source: Luna Innovations

Luna Innovations has leading experts in controlling the polarization and phase of light. The vast portfolio of polarization controllers can be credited to the established legacy and expertise of Luna’s two acquisitions: General Photonics (GP) and New Ridge Technologies (NRT). With Luna’s focus on fiber optic test and measurement, GP’s specialization in polarization control, and NRT’s legacy in both of these areas, the merging of these industry-titans allows Luna Innovations to lead in polarization control, management and measurement. Contact [Luna Innovations](#) today to learn more.

LUNA INNOVATIONS

301 1st Street, SW, Suite 200
Roanoke, VA 24011
Tel: (540) 769-8400

GLOBALSPEC

257 Fuller Road
Suite NFE 1100
Albany, NY 12203
Tel: (518) 880-0200

ABOUT LUNA INNOVATIONS

Luna Innovations Incorporated (NASDAQ: LUNA) is a leader in optical technology, committed to serving its customers with unique capabilities in high-performance, fiber-optic-based sensing, measurement, testing and control products for the aerospace, transportation, infrastructure, security, process control, communications, silicon photonics, defense, and automotive industries, among others. Enabling the future with fiber, Luna’s business model is designed to accelerate the process of bringing new and innovative technologies to market.